

## Setting up a Michelson interferometer on the laser optics base plate

### Objects of the experiment

- Assembling a Michelson interferometer
- Observing the interference pattern

### Principles

Interferometry is an extremely precise and sensitive measuring method for determining e.g. changes in lengths, layer densities, refractive indices and wavelengths. The Michelson interferometer belongs to the family of two-beam interferometers. It operates on the following principle:

The coherent light beam supplied by a suitable source is split into two parts by an optical component. These partial beams travel along different paths, are reflected into each other and channeled to another optical component, where they are combined and superimposed. The result is an interference pattern. If the path length of one of these partial beams, i.e. the product of the refractive index and geometric path, changes, this produces a phase shift with respect to the undisturbed beam. This in turn causes a change in the interference pattern, which allows us to draw conclusions about the changes in either the

refractive index or the optical path when the respective other quantity remains constant.

This means that, when the refractive index remains constant, we can determine differences in the geometric path, e.g. changes in the dimensions of materials due to heat or the effects of electric or magnetic fields. On the other hand, if the geometric path remains constant, we can determine refractive indices as well as quantities and influences which affect the refractive index such as changes in pressure, temperature or density.

The Michelson interferometer is suitable e.g. for demonstrating the effects of mechanical shocks and air streaking on the laser optics base plate. When setting up assemblies for making holograms, this arrangement can help you to recognize and eliminate disturbances.




Fig. 1: Photograph of the interference pattern on the translucent screen.

**Apparatus**

1 laser optics base plate . . . . .	473 40
1 He-Ne laser, linearly polarized . . . . .	471 840
1 laser support . . . . .	473 41
4 optics bases . . . . .	473 42
1 beam divider . . . . . e.g.	473 432
1 holder for beam divider . . . . .	473 43
2 planar mirrors with fine adjustment . . . . .	473 46
1 spherical lens, $f = 2.7$ mm . . . . .	473 47
1 translucent screen . . . . .	441 53
1 saddle base . . . . .	300 11
1 wooden ruler . . . . .	311 03

**Setup and Carrying out the Experiment**

*Note: optical components with damaged or dirty surfaces can cause disturbances in the interference pattern.*

*Handle the planar mirror, beam divider and spherical lens carefully, store them free of dust and do not touch them with your bare hands.*

Fig. 2 shows the setup of the Michelson interferometer on the laser optics base plate. To set up the experiment correctly, you must carry out the following steps:

**Laser optics base plate and laser:**

- Pump up the air cushion.
- Place the laser optics base plate **(a)** with air cushion horizontally on a sturdy laboratory bench.
- Mount the laser on the laser support and place it at the left edge of the base plate.
- Connect the laser and switch it on.
- Loosen the three lock nuts of the adjusting screws on the laser support.
- Using the adjusting screws, adjust the height and inclination of the laser so that the beam travels perfectly horizontally about 75 mm above the base plate (there is still enough play for subsequent adjustment). Measure the spacing with the ruler.
- Tighten the lock nuts.

**Beam divider:**

*The reflected and transmitted partial beams should have similar intensities.*

*When using the variable beam divider (473 435), make sure that the laser beam strikes the beam divider more or less in the center.*

- First make sure that the beam divider **(b)** reflects the laser beam horizontally. To do this, place the beam divider with optics base in the beam path at the opposite end of the laser optics base plate and reflect the laser beam to a point next to the laser emission aperture.
- Correct the inclination of the beam divider, and thus the beam path, as necessary using the two adjusting screws on the rod.
- Finally, place the beam divider in the beam path at an angle of  $45^\circ$  as shown in Fig. 2. The partially transparent layer of the beam divider should face the laser.

**Planar mirrors:**

*Notes:*

*It is easier to adjust the setup in a somewhat darkened room.*

*In addition to the main beams, the multiple reflections also produce so-called parasitic partial beams of low intensity. These are subsequently screened out by the lens holders, and can thus be ignored in subsequent adjustment.*

*The quality of the laser beam is impaired when the partial beams reflected by the beam divider are reflected directly into the emission aperture of the laser.*

- Place planar mirror **(c)** so that the laser beam strikes it in the center.
- By turning the optics base on the laser optics base plate and manipulating the adjusting screws on the back, align the planar mirror so that the beam is virtually reflected into itself and, after transmission through the beam divider, is incident at a point just above the emission aperture of the laser.

**Safety note**

The He-Ne laser fulfills the German technical standard "Safety Requirements for Teaching and Training Equipment – Laser, DIN 58126, Part 6" for class 2 lasers. When the precautions described in the Instruction Sheet are observed, experimenting with the He-Ne laser is not dangerous.

- Never look directly into the direct or reflected laser beam.
- Do not exceed the glare limit (i.e. no observer should feel dazzled).




Fig. 2: Setup of the Michelson interferometer on the laser optics base plate, top view

- a laser optics base plate
- b beam divider
- c, d planar mirrors with fine adjustment
- e spherical lens
- f translucent screen

- Fasten the translucent screen (f) in the base and set it up behind the laser optics base plate as shown in Fig. 2 so that the laser beam strikes it in the center.
- Place planar mirror (d) in the partial beam reflected from the beam divider (b) as shown in Fig. 2; position it about as far away from the beam divider as planar mirror (c).
- By turning the optics base on the laser optics base plate and adjusting the screws, align the planar mirror so that this partial beam is also virtually reflected into itself and is recombined with the first partial beam after transmission through the beam divider.
- Adjust the planar mirrors (c) and (d) using the adjusting screws so that the most intensive beams of the two reflection groups are completely coincident on the screen.

**Spherical lens:**

- Place the spherical lens (e) on the laser optics base plate between the beam divider and the translucent screen to widen the beam (the small opening of the lens holder must face toward the beam divider).
- Adjust the height and lateral position of the spherical lens so that the two partial beams pass through it axially.

**Fine adjustment:**

If you do not yet see a pattern of lines on the translucent screen:

- Change the beam path by slightly changing the alignment of the beam divider or the planar mirrors; readjust the spherical lens as necessary.

The more the partial beams run in parallel between the beam divider and the screen, the wider and farther apart the interference lines are.

- Adjust the interference pattern so that it is easy to observe by slightly changing the alignment of the beam divider or the planar mirrors.

If you cannot achieve a satisfactory image by fine adjustment, repeat the interferometer adjustment procedure from the beginning.

*The interference pattern is much brighter and easier to observe when the laser is switched to an output power of 1 mW. As this can change the beam path slightly, you may need to adjust the beam path or the position of the spherical lens.*

**Measurement example**

Fig. 1 on the title page shows a photograph of the interference pattern on the translucent screen.