Physik

Chemie · Biologie

Technik

LEYBOLD DIDACTIC GMBH

12/88 -Sf-

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The Franck-Hertz tube is used to demonstrate the discontinuous emission of energy of free electrons to mercury atoms (electron collision experiment by Franck and Hertz).

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Note:

You can find information on using the Franck-Hertz tube (555 85) with the Franck-Hertz supply unit (555 88), which supplies all necessary operating voltages in a form usable for conducting experiments, in the corresponding Instruction Sheet (point-by-point recording of the Franck-Hertz curve, oscilloscope display, recording with an XY-recorder and computer-assisted recording).

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1 Safety Precautions

- Vacuum tube with mercury filling (approx. 5 g). Danger of implosion by impact, fall or similar. In case of breakage, collect and eliminate the mercury (using mercury adsorbens, 306 83).
- Do not apply any voltages to the cold tube (risk of short-circuit by metallic mercury between the electrodes).
- If gas discharge occurs, as indicated by abrupt current rise, immediately set both grid voltages to 0.
- Handle the hot pieces of equipment using a heat-resistant cloth or gloves.
- Ensure via control of the oven temperature (i.e. the oven voltage) that the maximum temperature limit for the Franck-Hertz tube is not exceeded (permanent operation 200 °C, short-term peak value 220 °C).

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2 Description, Technical Data

2.1 Franck-Hertz tube (555 85)

- ① Vacuum tube with cylindrically arranged electrode system, discharge pins for leakage currents and with approx. 5 g of metallic mercury;

Heater voltage and current: 3.15 V, 0.4 A
(adapted to a 6.3 V power supply by a resistor in socket 555 86)

Control voltage U_1 between cathode and grid g_1 : 0 to 4 V d.c.

Acceleration voltage U_2 between grid g_1 and grid g_2 : 0 to 40 V d.c.

Countervoltage U_3 between grid g_1 and collector A: 1.5 V approx.

Max. admiss. temperature: 200 °C (cont. duty)
220 °C (short-term)

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
2.2 Socket for Franck-Hertz tube (555 86)

- ② Socket with 6 colour-coded, marked connecting leads in shielded 5-core cable and separate shielding for the connecting lead to the collector A; with built-in dropping resistor to reduce a supply voltage of 6.3 V to the filament voltage of 3.15 V required for the tube.

- ③ Copper tube (100 mm long, 35 mm dia.) to adapt the Franck-Hertz tube to the diameter of the electric oven (555 81/82) and for screening against stray electric fields by earthing them via a copper strand with 4-mm plug. with blind hole (3.1) for the measuring point of a temperature probe (e.g. 383 02)

Fig. 2

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or (see Fig. 4.3)

AC mains voltage (only 220 V AC) via
Safety connecting box 502 06
connected to digital temperature
controller and indicator 666 198
with temperature probe 383 02

Recommended for temperature control (when working without temperature regulation)

Temperature probe 383 02
Digital thermometer 666 190

Voltage sources for 0 to 30 V DC, 6 V AC, 1 V to 1,5 V DC, 1,5 V DC:

either (see Fig. 3.1)

1 Stabilized power supply unit, 300 V AC .. 522 35
with voltage divider for 0 to 30 V DC, e.g.
Measuring resistor, 10 k Ω 536 19
Plug-in potentiometer, 4,7 k Ω 577 95
or
Plug-in motor potentiometer, 4,7 k Ω 581 49
with DC voltage source, approx. 1,5 V

1 Battery, 1,5 V from 503 11
1 Battery holder 576 86

or (see Fig. 3.2)

1 Stabilized power supply unit, ± 15 V DC .. 522 30
1 Transformer 6 V, 12 V; 30 VA 562 73
2 Batteries, 1,5 V from 503 11
2 Battery holders 576 86
1 Plug-in potentiometer, 220 Ω 577 90

Ammeter for 10⁻⁹ A DC,

either (see Fig. 3.1)

1 Current measuring amplifier D 532 00
with analogue voltmeter,
measuring range 3 V DC e.g. 531 781
and adapter BNC/4 mm, single-pole 501 09

or (see Fig. 3.2)

1 Electrometer amplifier 532 14
with resistor 1 G Ω 577 02
and voltage source 12 V AC e.g. 562 73
and analogue voltmeter,
measuring range 3 V DC e.g. 531 781

2 Voltmeters for 3 V AC and 30 V AC, e.g.
Multimeter (for U₁) 531 55
Demo-multimeter (for U₂) 531 91

1 Rastered socket panel A4 576 74

1 Set of 10 bridging plugs 501 48

For recording the Franck-Hertz curve

1 XY recorder e.g. 575 662

Recommended in the event that I_A is applied via the I measuring amplifier D to the Y input of the recorder (see Fig. 3.1):

1 Smoothing capacitor 2200 μ F 538 31

Recommended, if the accelerating voltage U₂ is not adjusted via the motor-potentiometer (581 49) but

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manually via potentiometer (577 95):

- 1 Smoothing capacitor, parallel to recorder input X, e.g. plug-in capacitor, 100 μF 578 39

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3.2 Preparing the experiment assembly for recording the Franck-Hertz curve

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! Important for experiment setup in accordance with Fig. 3.1:

- Since the stabilized power pack, 300 V DC (522 35) supplies voltage dangerous to contact, only carry out work on the circuit when the instrument has been switched off.
- Do not apply voltages higher than 30 DC to the rastered socket panel.


Fig. 3.1

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Fig. 3.1/3.2 Assemblies for measuring the Franck-Hertz curve using a recorder; For a point-by-point measurement do not connect the recorder as indicated by the dotted lines.


Fig. 3.2

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Fig. 4
Cylinder oven voltage supply for heating the Franck-Hertz tube to an operating temperature of approx. 190 °C.


Fig. 4.1

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Fig. 4.1/4.2
Temperature setting via oven operating voltage; temperature measurement by a temperature probe protruding through the oven hole into the blind hole (3.1) of the copper cylinder.


Fig. 4.2


Fig. 4.3

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Fig. 4.3
Temperature control via a temperature probe and temperature measuring and control unit (666 198).

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Set up the assembly according to fig. 3 and 4; do not connect the recorder yet.

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Carry out the following steps in the indicated sequence before connecting the recorder:

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- a) Operate the oven (with copper cylinder and tube inserted) at first for 5 min at 220 V AC (or 110 V) and then for 10 min at 90 V AC (or 45 V) or via the temperature measuring and control instrument for 20 min at a set point of 190 °C.

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Fig. 5
 I_A as a function of the accelerating voltage U_2
 (recorder diagram).

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- b) Switch on the cathode supply voltage.
- c) Only necessary when stabilized power pack (522 35) according to Fig. 3.1 is used:
 Turn the potentiometer (a) fully clockwise; apply a voltage U_2 of e.g. 30 V via potentiometer (c) and immediately reset potentiometer (a) to 0.
- d) After approx. 1 minute's cathode heating, slowly increase voltage U_2 via potentiometer (a); if the current rises slowly and continuously or it shows the typical Franck-Hertz characteristic, go on with step g).
- e) In case of abrupt current rise (gas discharge due to cold tube) immediately reset U_2 to 0. Increase oven voltage by 5 V or increase the set point by approx. 2 °C; after 5 min proceed as under section d).
 If necessary, keep on increasing the oven temperature until current I_A shows the characteristic described under section d).
- f) If the current I_A increases only very slowly when running through voltage U_2 , reduce the oven voltage by 5 V or reduce the set point by 2 °C. After cooling down the unit for 10 min, proceed as under d). If necessary, keep on reducing the temperature until I_A shows the characteristic described under section d).
- f) Run through U_2 at different control voltages (increase U_1 by steps of 0.1 V via potentiometer (b)) until the current shows at least 4 distinct peaks.
- h) If necessary, connect the recorder;
 calibration:
 X-axis: 1 V cm^{-1} ; var.
 Y-axis: 1 V cm^{-1} ; var.

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After finishing the experiment, at first switch off the operating voltage and then remove the tube from the oven.