Physics

Chemistry · Biology


Technology



Lehr- und Didaktiksysteme LD Didactic GmbH

Leyboldstrasse 1 · D-50354 Huerth

06/05-W97-Sel



Instruction sheet 562 11

U core with yoke (562 11) Clamping device (562 12)

1 U core (1a) with yoke (1b), fixing bolt (1c)
2 Clamping device

1 Description

The U core with yoke (562 11) and the clamping device (562 12) are intended for the assembly of the demountable transformer with the coils indicated in chapter 3.

2 Technical data

Cross section: $4 \text{ cm} \times 4 \text{ cm}$


Width of U core and

length of yoke: 15 cm Heigth of U core: 13 cm

Material: Iron, laminated

Max. relative permeability μ_{r}

at the initial curve: approx. 2200 at $H = 120 \text{ A m}^{-1}$



Hysteresis curve of the U core with yoke

3 Suitable coils

Low voltage coil <i>N</i> = 250	562 13
Low voltage coil N = 480, I = 10 A	562 131
Low voltage coil N = 500	562 14
Low voltage coil <i>N</i> = 1000	562 15
High voltage coil N = 10000	562 16
High voltage coil N = 23000	562 17
Extra-low voltage coil N = 50	562 18
Mains coil 230 V	562 21
Mains coil 115 V	562 22

4 Use

- Mount primary and secondary coil.
- Put unpainted side of the yoke on the U core.
- Clamp yoke using the clamping device.

