

Physik


Chemie · Biologie

Technik

LEYBOLD

LEYBOLD DIDACTIC GMBH

9/85 — Kp —

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Coil with 250 Turns  
 Coil with 500 Turns  
 Coil with 1000 Turns  
 Coil with 10000 Turns  
 Extra-low Voltage Coil  
 with 50 Turns

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### 1 Safety Notes

- Do not exceed the maximum permissible current for continuous and short-term operation!
- After an overload during short-term operation, allow the coils to first cool down before switching on the current again!
- If the coils are used in low-voltage and high-voltage circuits, we recommend that laboratory safety cables (500 60 — 84) be used; if necessary, refit equipment without safety sockets with safety adapter sockets (500 95 — 98).
- Connect the primary coil of a transformer via safety connecting box (502 06) to the AC voltage supply; make alterations to the primary circuit of a mains transformer only with the primary voltage switched off.
- Before making alterations to the secondary circuit of a split low/high-voltage transformer, switch off the primary voltage.
- Do not obstruct the ventilation slits in the coil casing.

### 2 Technical Data

Coil Cat. No.	Number of turns	Permanent load current A	Resistance Ω	Inductance without iron core	With thickness mm dia.
562 13	250	5	approx. 0.6	0.0022	1.5
	Primary coil for mains transformer 110 V and for electromagnet				
562 14	500	2.5	approx. 2.5	0.009	1.0
	Primary coil for mains transformer 220 V				
562 15	1000	1.25	approx. 9.5	0.036	0.7
562 16	10000	0.1	approx. 1350	4.5	0.2
562 18	approx. 50	10	approx. 0.08	0.0001	2.0
	Secondary coil for low-voltage transformer with 562 14 or 13 as primary coil on U-core with yoke and clamping device, 562 11/12				

Coil casing of plastic material (approx. 8 cm x 8.5 cm x 11 cm) with ventilation slits;

direction of winding indicated on top and bottom of enclosure;

beginning, ending and centre tap accessible via 4 mm safety sockets.