11/94-Sf-

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The apparatus is used to demonstrate the Hall voltage U_H on a silver or tungsten strip bearing a current I in a magnetic field B acting at right angles to the direction of current.

When using the Hall effect apparatus, silver (586 81), quantitative investigation of the relationships

$$U_H = f(I) \text{ and } U_H = f(B)$$

via the proportionalities

$$U_H \sim I \text{ and } U_H \sim B$$

gives

$$U_H = \text{const.} \cdot I \cdot B$$

The theoretically derived formula for the Hall voltage of a strip-shaped conductor (of thickness d) made of material with a charge carrier concentration n is thus confirmed:

$$U_H = \frac{1}{n \cdot e} \cdot \frac{1}{d} \cdot I \cdot B \quad (I)$$

The material-dependent factor $\frac{1}{n \cdot e}$ is designated as the Hall constant R_H .

In equation (I), all quantities except n are accessible to measurement) so that the Hall effect enables the charge carrier concentration to be determined by experiment.

The direction of the Hall voltage in silver indicates negative charge carriers. The result is in agreement with the concepts of the model of the free electron gas. According to this, the most weakly bonded electrons (valency electrons) - for silver e.g. 1 valency electron per atom - are moving freely within the metal.

The limits of this model are shown by the so-called "abnormal Hall effect" of tungsten. Experiments carried out with the Hall effect apparatus, tungsten (586 84) under identical conditions give the following result: The Hall voltage in tungsten has the same magnitude but the opposite direction as in silver. This can be explained by the "energy band diagram". So-called hole electrons or holes are assigned to empty positions near the upper edge of an otherwise filled energy band. As can also be shown theoretically, such holes behave relative to electric or magnetic fields in such a way that they would seem to have positive charges.

Literature:

Physics Experiments, Volume 3 (599 942)

New Physics Leaflets for Colleges and Universities, Volume 1 (599 952)

New Physics Leaflets for Colleges and Universities, Volume 2 (599 892)

^{*}) The elementary charge e can e. g. be determined using the Millikan apparatus (559 41/42), unless it is given as natural constant.


Fig. 1

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1 Safety notes

For transverse currents over 15 A or magnet currents above 5 A, only switch on the device briefly (overheating of leads or overloading of the coils, which are designed for a maximum load of 5 A).

In the transverse current circuit, use cables which are rated for a maximum load of 20 A (e.g. connecting leads 501 20 ff or safety connecting leads 500 610).

Protect the experiment setup from drafts while measuring the Hall voltage.

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2 Description and technical data

- ① Strip type conductor of silver (included in 586 81) or tungsten (included in 586 84);
strip thickness $d = 5 \cdot 10^{-5}$ m
- ② Sockets for transverse current I , max. admissible intensity of current I : 22 A DC (for short periods).
- ③ Pair of sockets for tapping the Hall voltage U_H ; with polarity designation (for indicating a positive Hall voltage).
Magnitude of Hall voltage U_H : 10^{-6} V.
- ④ Adjusting knob for built-in 5 Ohm potentiometer for zero point adjustment
- ⑤ Stand rod to hold the apparatus in the bore hole of the U-core (from 562 11), to which the electromagnet for generation of the homogeneous magnetic field is fitted.
Required field strength B : 0.1 T to 0.9 T

Dimensions of panel: approx. 13 cm x 16 cm x 0.2 cm

Weight: approx. 0.4 kg

3 Operation

3.1 Circuit assembly; equipment and measuring instruments

For both Hall apparatus the circuit is assembled according to Fig. 2, using the following equipment and measuring instruments for the transverse current I , the magnetic field and the Hall voltage.

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Fig. 2 Experiment assembly for Hall Effect

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or

Regulated power supply unit, 0 to 12 V, 20 A 522 47

Current change via two rheostats 11Ω (537 26) connected in parallel or via a piece of constantan wire, 0.5 mm dia (550 44).

Ammeter, range 30 A DC, e. g.

Demo-multimeter 531 911

Shunt resistor up to 30 A 531 92

3.1.2 Magnetic field

Electromagnet, assembled from

U-core 562 11

Pair of bored pole pieces 560 31

2 coils, 250 turns 562 13

DC voltage source, adjustable up to approx. 20 V, 10 A,

Variable low-voltage power supply 522 39

Ammeter for magnet current, range 6 A DC or 10 A DC, e. g.

AV-Meter 531 94

Measuring device for magnetic field strength B , e. g.

Tangential B -probe 516 60

Teslameter 516 62

3.1.3 Hall voltage U_H

Voltage-sensitive measuring instrument, range 10^{-5} V to 10^{-6} V DC, e. g.

Microvoltmeter 532 13

Recommended:

Analog measured-value indicator, e.g.

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
4 Notes on experiments

Fit the Hall effect apparatus according to Fig. 2 into the electromagnet whose pole pieces are placed right against the panel so as to keep the air gap where the silver or tungsten is positioned as narrow as possible.

Important: Before recording the $I_M - B$ calibration curve and prior to every determination of B from this curve, demagnetize the iron of the electromagnet (allow a current of approx. 5 A AC to flow through the coils for a short time, then steadily reduce it to zero).

Determine the calibration curve $B = f(I_M)$ without the Hall effect apparatus using exactly the same pole-piece spacing as is required by the panel in experiments on the Hall effect.

Put the *initially unconnected, voltage-sensitive measuring instrument for the Hall voltage U_H* (e.g. microvoltmeter, 532 13) into operation according to the respective instruction sheet and set the display to zero with Reset.

Before exposing the Hall effect apparatus to the magnetic field, adjust the zero point: Apply a transverse current I of e.g. 10 A and set the indicator of the meter for measuring the Hall voltage U_H to zero using the adjusting knob ④. If the display changes after switching off, switch the transverse current back on and repeat the zero-point adjustment.

To demonstrate the normal and the abnormal Hall effect (verification of magnitude and direction of Hall voltage), carry out experiments using the Hall effect apparatus (silver) and the Hall effect apparatus (tungsten), under identical conditions, with max. operating data (magnet current $I_M = 10$ A, transverse current $I = 20$ A).

For demonstrating the proportionalities $U_H \sim I$ and $U_H \sim B$ and for precise determination of the Hall voltage U_H , the Hall effect apparatus (silver) (586 81) is most suitable.

Quantitative experiments with the tungsten apparatus call for special care and dexterity on the part of the experimenter.

With switched-on transverse current, air circulation may cause considerable zero-point fluctuations (thermal e.m.f. on the measuring contacts for Hall voltage).

Due to the higher electric resistance of tungsten, the thermal effects and hence the zero-point fluctuations are higher than with silver.

Carry out at least five measurements to determine the Hall voltage.

5 Example of measurement using the Hall effect apparatus, silver, 586 81 (see Fig. 3)

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Tabellenwert: $R_H = 8,9 \cdot 10^{-11} \text{ m}^3 \text{ C}^{-1}$

Hall voltage U_H as a function of the magnetic flux density B at a transverse current I of 20 A; the slope $\Delta U_H/\Delta B$ is given for calculating the Hall constant $R_H = \frac{1}{n \cdot e} \cdot \frac{1}{d} \cdot I \cdot B$

$$R_H = \frac{\Delta U_H}{\Delta B} \cdot \frac{d}{I} = \frac{29,4 \cdot 10^{-6} \text{ V T}^{-1} \cdot 5 \cdot 10^{-5} \text{ m}}{20 \text{ A}} = 7,35 \cdot 10^{-11} \text{ m}^3 \text{ C}^{-1}$$

Table value: $R_H = 8,9 \cdot 10^{-11} \text{ m}^3 \text{ C}^{-1}$