
Physik

Chemie · Biologie

Technik

LEYBOLD DIDACTIC GMBH

9/95-Wit-

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Apparatus Collection for Laser Optics (Interferometry and Holography)

Laser optics base plate (473 40)
Laser mount (473 41)
Optics base (473 43)
Holder for beam divider (473 43)
Beam divider 50% (473 432)
Variable beam divider (473 435)
Planar mirror with fine adjustment (473 46)
Spherical lens, $f = 2.7$ mm (473 47)
Film holder (473 44)
Object holder (473 45)
Fine-adjustment mechanism (473 48)
Evacuatable chamber (473 485)

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This apparatus collection allows you to assemble a variety of interferometer types and to conduct experiments with these arrangements. In addition, this collection can also be used to create and reconstruct several different kinds of holograms.

The He-Ne laser, linearly polarized (471 840), is required as a light source.

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1 Safety notes

- The He-Ne laser, polarized (471 840) fulfills the German technical standard "Safety Requirements for Teaching and Training Equipment - Laser, DIN 58126, Part 6" for class 2 lasers. When the precautions described in the respective Instruction Sheet are observed, experimenting with the He-Ne laser is not dangerous.
- Never look directly into the direct or reflected laser beam.
- Do not exceed the glare limit (i.e. no observer should feel dazzled).
- The laser optics base plate is heavy (30 kg)! Take this into consideration when transporting this device.
- Place the laser optics base plate on the air cushion and protect it from mechanical bumps or shocks during the experiment.
- Store the spherical lens, beam dividers and planar mirror protected from dust; do not touch the surfaces.
- Handle the beam dividers only by the matt end faces at all times.
- Store the planar mirror with its protective cap.
- Carefully clean dirty optical components before the experiment; use e.g. lens tissue (305 00). Alcohol is a suitable cleaning fluid (see section 4.1/4.2).


Fig. 1
Laseroptik-Grundplatte
Laser optics base plate

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
Fig. 2
Laserträger
Laser mount

2 Description, technical data, scope of supply

2.1 Laser optics base plate (473 40)

- ① Laser optics base plate, granite, vibration-damped, with hinge (1.1) for attaching cover ③
Dimensions: 60 cm x 30 cm x 6 cm
Weight: approx. 30 kg
- ② Air cushion (inflatable rubber ring) for vibration isolation
Diameter (inflated): approx. 65 cm
- ③ Cover for minimizing convection when creating holograms, with hinge (1.1) for attachment to the laser optics base plate.
Dimensions: approx. 14 cm x 38.5 cm x 28 cm (H x W x D)

2.2 Laser mount (473 41)

For mounting the He-Ne laser, linearly polarized (471 840). Height and inclination adjustable, stable three-point stand.

- ④ Front mounting plate
- ⑤ Rear mounting plate
- ⑥ Knurled screws with lock nuts for height and inclination adjustment of laser;
Height adjustment range of laser beam approx. 71 mm to 79 mm
- ⑦ Centering taps for correct positioning
- ⑧ Square nuts to fit mounting groove of laser
- ⑨ Allan screw (size 5 key) for fixing in mounting groove

Not shown: Allan key (size 5 key)

Diameter of mounting plates: 8 cm each
Total weight: approx. 1.36 kg


Fig. 3
Optik-Fuß
Optics base

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2.3 Optics base (473 42)

- ⑩ Base plate with round recess for reducing the distance between midpoints when using multiple optics bases; with stable three-point stand
- ⑪ Clamping column for rods of up to 14 mm dia.
- ⑫ Clamping screw
Height: 4 cm
Outside diameter of base plate: 4 cm
Weight: approx. 310 g

2.4 Beam divider 50 % (473 432)

For dividing the light beam into a reflected and a transmitted component beam;
for experiments on interferometry

Semitransparent glass pane with beam divider coating on one surface
Intensity ratio of reflection to transmission at an angle of incidence of 45 ° approx. 1 : 1
For mounting in holder for beam divider (473 43)
Dimensions: 10 mm x 30 mm x 6 mm (H x W x D)

2.5 Variable beam divider (473 435)

For dividing the light beam into a reflected and a transmitted component beam;
for experiments on holography and interferometry

Semitransparent glass pane with a beam divider layer which varies continuously across the surface of the glass
Intensity ratio of reflection to transmission variable over the width
For mounting in holder for beam divider (473 43)
Dimensions: 10 mm x 58 mm x 6 mm (H x W x D)

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
Fig. 7
 Planspiegel
 Planar mirror

2.6 Holder for beam divider (473 43)

- ⑬ Support with groove for accommodating beam dividers (473 432/435), vertical adjustment of beam dividers possible with grub screws ⑯
- ⑭ Clamping frame for holding beam dividers
- ⑮ Knurled screws for easy insertion and removal of beam dividers
- ⑯ Grub screws with Allan socket (size 2 key) for vertical adjustment of support ⑬; see section 3.2.2
- ⑰ Rod (dia. 12 mm) for mounting in optics base (473 42)

Not shown: Allan key (size 2 key)

2.7 Planar mirror with fine adjustment (473 46)

- ⑱ Planar mirror, dia. 4 cm
 Mirror surface: aluminum with protective SiO₂ coating;
 adjustable via knurled screws ㉔
- ⑲ Plastic holder (dia. 5.3 cm)
- ㉔ Knurled screws for fine adjustment of planar mirror in the optical arrangement
 (top knurled screw for vertical adjustment
 bottom knurled screw for horizontal adjustment)
- ㉕ Grub screw (size 2 key) for pretensioning of the adjustment mechanism of the planar mirror (do not change; see section 4.1, Maintenance)
- ㉖ Rod (dia. 12 mm) for mounting in optics base (473 42) or in fine-adjustment mechanism (473 48)
- ㉗ Protective cap for planar mirror

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
Fig. 9
Filmhalter
Film holder

2.8 Spherical lens, $f = 2.7 \text{ mm}$ (473 47)

- ②④ Lens body, with antireflective coating
Focal length: 2.7 mm
Diameter: 5 mm
- ②⑤ Plastic holder, dia. 2.2 cm
Beam opening dia. 3 mm
- ②⑥ Locking bolt with integrated beam opening for securing lens body, fixed with grub screws (26.1) (size 2 key);
Bolt with grip notch for easy removal from holder ②⑤ (required for cleaning the lens body, see section 4.2, Maintenance)
Beam opening dia. 2 mm
- ②⑦ Rod (dia. 12 mm) for mounting in optics base (473 42)

Not shown: Allan key (size 2 key)

2.9 Film holder (473 44)

- ②⑧ Aluminum film holder frame, for mounting flat film cut to size 42 mm x 51 mm
Light aperture 37 mm x 47 mm (H x B)
- ②⑨ Groove for fingernail (to aid correct insertion of film in darkened room)
- ③⑩ Clamp for securing film
- ③① Knurled screw for opening and closing clamps ③⑩
to open, tighten the knurled screw
to close, loosen the knurled screw
- ③② V-grooves for positioning pins ③⑤ for reproducible and stable positioning of the film holder frame in the experiment assembly
- ③③ Magnet for attaching the film holder frame ②⑧ to the film holder support ③④ (contact-free magnetic holder)
- ③④ Steel film holder support, for accommodating the film holder frame ②⑧
- ③⑤ Positioning pins for supporting the film holder frame in the V-grooves ③②
- ③⑥ Rod (dia. 12 mm) for mounting in optics base (473 42)


Fig. 10
Objekthalter
Object holder

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2.10 Object holder (473 45)

- ③⑦ Steel object holder plate with stage, with 2 threaded holes (37.1), (37.2) for mounting post ③⑧ at different positions
Plate thickness: 10 mm
Plate area: 57 mm x 50 mm
Stage: 23 mm x 12 mm
- ③⑧ Post with knurled knob for easy repositioning; with holder arm (38.1) for fixing objects
Holder arm with rubberized bottom edge to prevent objects from slipping
Clamping height: approx. 5 cm
- ③⑨ Rod (dia. 12 mm) for mounting in optics base (473 42)

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2.11 Fine-adjustment mechanism (473 48)

For determining the wavelength of laser light by means of defined displacement of the planar mirrors (473 46 / 471 87) in the Michelson interferometer;
suitable for use on the laser optics base plate (473 40) and the interferometer base plate (471 85)
④① Reducing gear 100 : 1



Fig. 11
Feinstelltrieb
Fine-adjustment mechanism

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- (40.1) Double universal-joint coupling for connection to the joint head (41.1)
- (40.2) Adhesive magnetic strip for attaching gear unit to gear base ④②
- (40.3) Knob with crank
- ④① Displacement assembly with fixed micrometer screw
 - (41.1) Joint head for connection to double universal-joint coupling (40.1)
 - (41.2) Micrometer screw:
Advance: 0.5 mm / turn
Length of stroke: max 25 mm
Division: 0.01 mm
 - (41.3) Rider for horizontal and vertical attachment of planar mirror
 - (41.3.1) Horizontal mount for planar mirror with fine adjustment (473 46), for use on the laser optics base plate (473 40)
 - (41.3.2) Vertical mount for front silvered plane mirror (471 87), for use on the interferometer base plate (471 85)
 - (41.4) Knurled screw for clamping the elements in (41.3)
 - (41.5) Pull-back spring for rider (41.3)
 - (41.6) Rod (dia. 12 mm) for attachment in stand material, e.g. optics base (473 42) or magnetic lens base (471 86)
- ④② Gear base for reducing gear ④①
 - (42.1) Surface (11 cm x 3 cm) for mounting the reducing gear
(attached using adhesive magnetic strip (40.2))
 - (42.2) Rod (dia. 10 mm) for attachment in stand material, e.g. optics base (473 42) or magnetic lens base (471 86)

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For determining the index of refraction of gases using interferometric means

- ④③ Aluminum jacket
- ④④ Glass pane glued vacuum-tight
- ④⑤ Hose connectors (outside dia. 12 mm, inside dia. 6 mm) for connecting the vacuum pump (or gas supply)
- ④⑥ 2 rubber stoppers (4/8 mm) for sealing the hose connectors
- ④⑦ Rod (dia. 10 mm, length 3.6 cm) for mounting in optics base (473 42)

Internal chamber dimensions: length 50 mm, dia. 40 mm
Outside dimensions: length 58 mm, dia. 59 mm
Weight: approx. 260 g


Fig. 13
Befestigung der Abdeckhaube an der Laseroptik-Grundplatte
Attaching the cover to the laser optics base plate

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3 Operation

For detailed information on setting up and carrying out experiments, refer to the corresponding experiment descriptions.

3.1 Experiment preparation

3.1.1 Laser optics base plate (473 40)

Pump up the air cushion ② and place it together with the laser optics base plate ① in a horizontal position on a sturdy laboratory bench insensitive to shocks. During the experiment, protect the base plate from mechanical shocks and avoid air streaking in the setup, e.g. through breathing or drafts.

For assemblies for creating holograms: attach the cover ③ to the laser optics base plate as shown in Fig. 13 to protect the experiment zone from convection during exposure.

When attaching the cover for the first time, you may need to adjust the seat on the laser optics base plate. To do this, close the cover and loosen the mounting screws of the hinge on the rear of the cover slightly, so that the cover lies flush on the laser optics base plate. Then retighten the mounting screws.

3.1.2 Mounting and adjusting the He-Ne laser, linearly polarized (471 840) on the laser mount (473 41)

Attach the laser mount to the laser. To do this, slide the square nuts ⑧ into the mounting groove on the bottom of the laser, and attach the mounting plates ④ and ⑤ to the laser using the Allan screws ⑨ as shown in Fig. 14. Attach the front mounting plate ④ beside the laser emission aperture, and the rear mounting plate ⑤ as close as possible to the other end of the laser.

Place the laser with laser mount on the laser optics base plate as required for the experiment. Loosen the lock nuts of the three adjusting screws ⑥ of the laser mount. Using these adjusting screws, alter the height and inclination of the laser so that its beam is parallel to and about 7.5 cm above the base plate. Then retighten the lock nuts.

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3.2 Beam dividers

3.2.1 Mounting the beam dividers (473 432/435) in the holder (473 43)

Important: in order to prevent damage to the beam-divider coating, handle the beam dividers only by the matt end surfaces. Carefully clean dirty or dusty glass surfaces using lens tissue (305 00). Do not touch these surfaces with your bare hands.

Remove the clamping frame ⑭ from the support ⑬ by loosening the knurled screws ⑮. Place the beam divider in the center of the support, then reattach the clamping frame to the base plate.

3.2.2 Adjusting the inclination

Note: the inclination adjustment allows you to make sure that the partial beam reflected by the beam divider is parallel to the surface of the laser optics base plate. This step must always be carried out after placing the beam divider in the holder. However, we recommend checking the adjustment before every experiment, and correcting as necessary.

Mount the beam divider with holder in the optics base and place it on the laser optics base plate opposite the laser after its horizontal beam path has been set. Align the beam divider so that it reflects the laser beam to a point next to the laser emission aperture. Using the grub screws ⑯ **with Allan socket (size 2 key),adjust the inclination of the beam divider as follows:**

- If the reflection is *above* the laser emission aperture, first loosen the grub screw facing *away* from the laser somewhat, and then carefully tighten the grub screw facing *toward* the laser until the reflection is on the level of the laser beam emission opening.
- If the reflection is *below* the beam emission aperture, loosen the grub screw facing *toward* the laser and carefully tighten the grub screw facing *away* from it.

Important: the adjusted inclination is not fixed until both grub screws ⑯ have been tightened slightly.


Fig. 16
Schneiden des Holografie-Films
Cutting the holography film

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3.3 Mounting the film in the film holder (473 44)

Please refer to the corresponding experiments for more detailed information on setting up the experiment and creating holograms.

3.3.1 Film and chemicals

Required for creating holograms:

1 Holography film ¹⁾	473 442
1 Darkroom accessories	473 446
1 Photographic chemicals ²⁾	473 444
For creating phase holograms: 671 891	
1 Iron(III) nitrate nonahydrate, 250 g	
1 Potassium bromide (KBr), 50 g	672 491

suitable:

¹⁾Agfa-Gevaert 2NFXQ HOLOFI 8E75 T3 HD NAH

²⁾BW-paper developer Agfa Neutol, fixer Tetenal Superfix

The bleach bath required for preparing phase holograms consists of

100 g iron(III) nitrate nonahydrate
30 g potassium bromide (KBr)
1 l distilled water

3.3.2 Cutting the film to size

As the film is extremely insensitive at a wavelength of 505 nm, you can use a weak, dark-green (or blue-green) darkroom lamp (e.g. from the darkroom accessories (473 446)) as an orientation aid. Recommendation: if the film is not to be used for several weeks, do not expose it to light.

Handle the film carefully, so as not to damage the coating. Note which side is coated! The film is marked by a notch in the edge. If this notch is at the top right or bottom left, the coated side is facing away from you.

This film material consists of coated plastic sheets (sheet film) 10.2 cm x 12.7 cm, which must be cut to a suitable format of 42 mm x 51 mm (tolerance 1 mm). In complete darkness, remove the desired number of sheets and pack up the remaining film light-tight. Under the appropriate darkroom light, mark the points at which to cut e.g. with a fine felt-tip pen (see Fig. 16) and cut to size with a scissors.

Store the pieces you have cut to size in an absolutely light-tight container (e.g. from the darkroom accessories (473 446)) (note which side is coated; if necessary, mark it by cutting off a corner) and use them within the following weeks.

3.3.3 Using the film holder

Set up the experiment according to the respective experiment description. To insert the film, hold the optics base of the film holder with your finger and carefully remove the film holder frame 29 from the setup without the film holder support 34. Open the


Fig. 18

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clamp ③ of the film holder by tightening the knurled screw ③. Take a piece of film cut to size from the light-tight box in a darkened room and slide it into the film holder so that the coated side is facing toward the object when the holder is returned to the setup. To make correct insertion easier, hold your thumb on the groove ② and slide the film as far as your fingernail (see Fig. 17). Close the clamp by loosening the knurled screw completely and clamp the film evenly.

Then, hold the optics base of the film holder with your finger and replace the film holder in the setup.

3.4 Fine-adjustment mechanism (473 48)

The fine-adjustment mechanism is designed for use on the laser optics base plate (473 40) and the interferometer base plate (471 85).

3.4.1 Fine-adjustment mechanism on the laser optics base plate (473 40)

Assemble the desired experiment setup and mount the displacement assembly ④ in the optics base (473 42). Align the recess of the optics base as shown in Fig. 18 so that the apparatus stands securely.

Fasten the planar mirror with fine adjustment (473 46) in the horizontal clamp (41.3.1). Mount the gear base ④ in the optics base and attach the reducing gear ④ to the gear base using the magnetic strip. Carefully link the double universal-joint coupling to the joint head of the micrometer screw of the fine-adjustment mechanism.

Adjust the position of the reducing gear so that the coupling rod is neither fully extended nor completely collapsed, as otherwise the measurement can be falsified later by displacement of the fine-adjustment mechanism.

Keep the angles of the individual segments of the universal joint coupling as small as possible (not greater than 45 °). Next, check the alignment of the optical component and readjust as necessary.

Turn the reducing gear slowly and evenly by lightly placing your finger on the crank of the gear knob and displace the planar mirror in this manner by 0.005 mm per revolution.

As the reducing gear with the joint coupling has a "dead" travel, it must be turned approx. six turns after changing direction before the element is displaced. Before starting the measurement, turn the gear at least one full turn using the knob. If the movement of the optical component is "jerky", lubricate the slide bush of the fine-adjustment mechanism.

3.4.2 Fine-adjustment mechanism on the interferometer base plate (471 85)

Mount the displacement assembly ④ and gear base ④ on one magnetic lens base (471 86) each. Attach the front-silvered plane mirror (471 87) in the vertical holder of the displacement assembly (see Fig. 19).

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3.5 Evacuatable chamber (473 485)

Seal one of the hose connections of the evacuable chamber tightly using a rubber stopper ④6. Then mount the evacuable chamber on an optics base and place it in one of the beam paths so that the partial beam passes through it axially. Make sure that no reflections of the beam from the glass surfaces strike the emission aperture of the laser, as this can effect the quality of the laser beam.

Additionally required for evacuating the chamber:

1 Hand vacuum and pressure pump	375 58
1 Saddle base, V-shape	300 02
1 Universal clamp S	666 555

Connect the vacuum pump to the other hose connector of the evacuable chamber as shown in F.g 20; use the correct hose adapter and plug it into the hose connector. Assemble the strain relief using the stand base and the universal clamp S and fasten the hose beside the laser optics base plate so that the measurement cannot be falsified by twisting or shifting the evacuable chamber.

4 Maintenance

4.1 Planar mirror with fine adjustment (473 46)

If the planar mirror cannot be satisfactorily adjusted using the knurled screws, reset the pretensioning of the adjustment mechanism using the grub screw ②1 with Allan socket (size 2 key) as shown in Fig. 21. The pretensioning is increased by screwing in the grub screw (turning clockwise).

4.2 Cleaning the spherical lens (473 47)

When creating holograms, it is important that the widened laser beams are as free as possible from high-contrast interference phenomena. These can be caused by dust particles, fingerprints or scratches in the surface of the lens. For this reason, you should store the spherical lens in a dust-free place and handle it carefully.

The lens body can be removed from the holder for cleaning. Be extremely careful when cleaning the lens.


Fig. 21 Vorspannungs-Einstellung der Justierklemme des Planspiegels
Setting the pretensioning of the adjusting mechanism of the planar mirror

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Important:

To prevent the little round lens body from becoming lost, you should carry out the following steps in or over a receptacle.
Be careful not to damage the surface of the lens. Do not touch it with your bare hands.
Use the lens tissue (305 00) for cleaning and handling the lens; alcohol can additionally be used to clean extremely dirty lenses.

4.2.1 Removing the lens

Place the holder 25 so that the locking bolt 26 faces upward. Loosen both grub screws (26.1) and remove the locking bolt. Hold the lens tissue by the opening of the holder and roll the lens into the lens tissue by tilting the holder.

4.2.2 Cleaning the lens

Clean the lens body with the lens tissue without touching it with your bare fingers. Also clean those parts of the holder and the locking bolt which are in direct contact with the lens using a lens tissue.

4.2.3 Inserting the lens

Place the lens body in the holder and insert the locking bolt in the holder (make sure it is the right way around). Then carefully tighten the grub screws.

To check how clean the beam is, place the spherical lens with optics base on the laser optics base plate in the beam path in front of the laser so that the small beam aperture faces the laser. Align the spherical lens vertically and horizontally so that the laser beam passes through it axially. In a darkened room, check the widened laser beam e.g. with a sheet of white paper. When the spherical lens is properly adjusted, the laser beam can show the following irregularities:

- A system of concentric rings

The lens body is provided with an antireflective coating; for manufacturing-related reasons, this has an uncoated band circling the lens body. This band is now in the optical axis, and the lens body must be turned in its holder somewhat.

Loosen the grub screws, pull out the locking bolt a bit, grasp the holder so that the locking bolt is facing up and carefully shake it. Then retighten the locking bolt and check the quality of the widened laser beam again.

- Point-type interference phenomena with irregular structure

These indicate that there is still dirt (e.g. dust) on the lens surface. In this case, the lens must be cleaned again. If the lens is only slightly dirty, it can still be used to create holograms and in interferometers, provided you are willing to accept the corresponding quality losses.

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